


## Great perspectives with microcrystals

Max IV workshop  
– structural biology

Oct. 21<sup>st</sup>, 2021



Poul Nissen

Danish Research Institute of Translational Neuroscience - DANDRITE  
Nordic-EMBL Partnership for Molecular Medicine & Lundbeck Foundation  
Center for Proteins in Memory – PROMEMO  
Danish National Research Foundation  
BRAINSTRUC  
Lundbeck Foundation  
Department of Molecular Biology and Genetics  
Aarhus University

DANDRITE



NORDIC EMBL  
PARTNERSHIP FOR  
MOLECULAR MEDICINE

EMBL



PROMEMO

Danmarks  
Grundforskningsfond  
Danish National  
Research Foundation



BRAINSTRUC



novobordisk  
fonden



INDEPENDENT RESEARCH  
FUND DENMARK



1

## MicroMAX, BioMAX and other biosample beamlines


### Unique resources:

MAX IV and ESS serve a general structural biology community using also e.g. cryo-EM possible interactions also with e.g. EMBL-Hamburg/DESY and European XFEL

### Structural biology 2021:

New, difficult structures will be based mostly on single-particle cryo-EM analysis  
Small globular structures are highly guided by AlphaFold – but ligand binding?

- Time-resolved cryo-EM is still in its infancy and based on cryo-conditions
- In situ cryo-electron tomography is limited to sub-micrometer vitrified samples/sections



### Large, important field for X-ray science:

Molecular and cellular structure/dynamics in biological context (time, temperature, 3D – living?)

Structure determination and time-resolved studies using serial crystallography/SAXS – *maintain crystallography/SAXS*

X-ray imaging of cells, organoids, native tissues – *support molecular cell biology, integrate with e.g. EMBL program*

Solution X-ray scattering studies is also an important complement to single-particle cryoEM

Integration of data and models – *support computing infrastructure for data analysis and modeling of complex dynamics*

Highly structured molecules/assemblies

Dynamic structures

Fully disordered molecules/assemblies



2

## Advantages of serial crystallography/microcrystals

- Less sensitive to mechanical stress and strain
- Fast diffusion rates in soaking
- Uniform activation of photocaged compounds
- Possibility for both room temperature and cryo-crystallography
- Manageable in small-scale format

### Potential applications

Time-resolved studies of e.g. photoactivated or diffusion-based reactions

Time-resolved studies of ligand binding reactions guiding drug design

(With microfocus: also X-ray absorption/spectroscopy/emission with subcellular resolution in imaging modalities)

### Challenges:

More diversity of photocaged compounds - chemistry

Integral technologies at beamlines for pump-probe/stopped-flow, laser activation, sample delivery

3

## Structural biology – connecting structure and dynamics to function


### SLC6 transporters

Neurotransmitter, monoamine  
and amino acid transporters  
E.g. serotonin, dopamine, glycine,  
GABA transporters, B<sup>0</sup>AT1  
Bacterial orthologs (amino acids)  
e.g. LeuT and MhsT

### Amino acid transporter MhsT

(*B. halodurans*)  
SLC6 family orthologue

- Occluded, inward-facing with bound Na<sup>+</sup> and Trp
- N-term bound or unbound/TM5 unwound or continuous
- Solvation and dynamics of Na<sup>+</sup> release coupled to TM5 and N-terminal dynamics



Malinauskaite *et al.* *Nature Struct Mol Biol.* 2014

Stolzenberg *et al.* *JBC* 2017

Gotfryd *et al.* *Nature Comm.* 2020

Focht, Neumann, Lyons *et al.* *EMBO J* 2021

4

## Article

# Structural insights into the inhibition of glycine reuptake




<https://doi.org/10.1038/s41586-021-03274-z>



Received: 29 May 2020

Accepted: 20 January 2021


Azadeh Shahsavar<sup>1,2</sup>, Peter Stohler<sup>2</sup>, Gleb Bourenkov<sup>2</sup>, Ivan Zimmermann<sup>4,5</sup>, Martin Siegrist<sup>2</sup>, Wolfgang Gube<sup>3</sup>, Emmanuel Pinard<sup>3</sup>, Steffen Sinning<sup>4</sup>, Markus A. Seeger<sup>4</sup>, Thomas R. Schneider<sup>2,5</sup>, Roger J. P. Dawson<sup>3,5,6</sup> & Poul Nissen<sup>1,5</sup>



Microcrystals  
in meso phase



Heatmap  
in MxCuBE





Shahsavar et al. Nature 2021



5

## Serial crystallography – synchrotron sources and XFEL

*B. halodurans* amino acid transporter MhsT-Met  
Inhibited state (non-substrate amino acid)  
Microcrystals, HiLiDe method




Serial Synchrotron X-ray  
crystallography (SSX):

- P14 beamline PETRA III  
- 5 x 10 μm<sup>2</sup>  
microfocus beam

- 1 μm sample displacement
- 0.2° oscillation
- 0.025 sec exposure time
- 100% transmission

Data collected from  
92 loops

Ellipsoidal




Trampari et al.  
in preparation

<b>Data collection</b>	P14
Space group	P2 <sub>1</sub>
Cell dimensions	
<i>a</i> , <i>b</i> , <i>c</i> (Å)	44.27, 215.32, 51.02
$\alpha$ , $\beta$ , $\gamma$ (°)	90, 90.111, 90
Resolution range (Å)	28.42 - 2.7/3.1/3.0 (2.838 - 2.74)
Total reflections	703132 (66641)
Unique reflections	24985 (245)
R-meas (%)	74.33 (953.8)
R-pim (%)	13.7 (184.1)
R-merge (%)	73.0 (935.4)
CC (1/2)	0.993 (0.354)
Mean I/ $\sigma$ I	5.50 (1.5)
Completeness (%)	81.2 (9.68)
Multiplicity	28.1 (26.3)
Wilson B	57.31





6


## Warming up for time-resolved crystallography



Samuel Hjort-Jensen  
(EU-ITN RAMP)



SERCA phosphorylation @RT  
[Ca<sub>2</sub>]E1.ADP-AlF<sub>4</sub>.2 Mg<sup>2+</sup>



Investing time in well-diffracting microcrystals

- Time-resolved binding/catalysis
- Drug discovery (intermediates)
- Need not be high resolution to restrain MD importantly